



# Spacer Use

Susan M. McKillop, MSNeD, BSN, RN, AE-C Asthma Educator

> Weill Cornell Medicine Pediatric Asthma Program




### Objectives

- lungs?
- Important things to know about using your spacer
- Using a spacer demo
  - Getting Ready
  - Breathing techniques
- Device Cleaning
- Refill/Replacing your inhaler

How does your asthma medication get into your



#### Where Does Your Asthma Medicine Go?



you can:



#### Things to know about using a spacer

- How to use the spacer with your asthma medications: - priming - breathing
  - cleaning
- Rinse the mouth/face after Inhaled corticosteroid (ICS) use
- Do not use a spacer with any dry powdered medications



#### How to Use a Spacer: <u>Getting Ready</u>

If you have not used the inhaler in a while, you may need to **prime** it

Take the cap off the inhaler and spacer.

Shake the inhaler hard 10 to 15 times before each use.

Attach the spacer to the inhaler.

Sit up straight and chin up

Hold the spacer away from your mouth and breathe out gently to empty your lungs.



### How to Use a Spacer: Breathe in Slowly

Put the spacer between your teeth and close your lips tightly around it.

If using face mask, secure mask by getting nose and mouth in mask space. Make sure fit is correct.

Are you still sitting up and chin up?

Spray one puff into the spacer by pressing down on the inhaler.

Breathe in slowly and breathe and deeply



### How to Use a Spacer: <u>Holding your Breath</u>

- For children that can not hold for a count of 10, have them breathe slow and regular for approximately 6 breaths.
- ▶ If you are using inhaled, quick-relief medicine (Albuterol, Ventilin, Pro-air), wait about 1 minute before you take your next puff.
- You do not need to wait a minute between puffs for your controller medicines.
- Put the caps back on the inhaler and spacer.

▶ Hold your breath as you count to 10, if you can. This lets the medicine reach deep into your lungs.

► After using your inhaler, if using an ICS, rinse your mouth with water, gargle, and spit. Do not swallow the water. This helps reduce the side effects from your medicine.



### Device Cleaning

To clean:

- For your inhaler:
  - Remove the metal canister from the L-shaped plastic mouthpiece •
  - Rinse <u>only</u> the mouthpiece and cap in warm water and dry
- For your spacer:
  - Cleaning should be done once a week.
  - Take the spacer apart and soak in warm soapy water.
  - Rinse and shake dry. (Never dry the inside of the spacer.
  - Let air dry on paper towel
  - A good way to keep both of these clean, place in a storage bag that can zip close



### Refill/Replacing your Inhaler

- Most inhalers come with counters on the canister. Keep an eye on the counter and replace the inhaler before you run out of medicine.
- Rule of thumb: if the canister reads "9", it is time to call the pharmacy for a refill!



## Contact Us

- Visit our website for links to more resources: <u>https://</u> <u>pediatrics.weill.cornell.edu/divisions-programs/general-academic-</u> <u>pediatrics/pediatric-asthma-program</u>

